

B302 Novel Binocular Biological Microscope

Product Description

B302 series biological microscope is a special type, it has beautiful drop-style appearance and perfect performance. B302 biology microscope adopts infinite optical system, and it is a professional labs microscope. The B302 series biological microscopy is cost-effective.

B302 series novel biological microscope includes 4 models (B302-1, B302-1TR, B302-2, B302-2TR). You can choose the biological microscope according to your demands.



Product Feature

Novel drop-style integrated appearance design

Standard binocular microscope and optional trinocular head

A pair of high eyepoint WF10X/20mm eyepiece

Low position coaxial coarse and fine focusing knobs provide comfortable long-term using

Standard cable wound device is convenient for using, storage

145*140mm double layer mechanical stage, X*Y moving 76*52mm



Product Parameter

Product Name	B302 Novel Binocular Biological Microscope	
Optical System	Infinite optical system	
Model	B302-1	B302-2
Observation Head	Binocular head ,inclined at 30°,360°rotatable. High-eyepoint WF10X/20mm eyepiece 48-75mm interpupillary distance	
Nosepiece	4 Holes	
Infinite Objective	EPlan 4X/10X/40X(S)/100X(S,Oil)	Plan 4X/10X/40X(S)/100X(S,Oil)
Platform	Size:145*140mm, travel 76*52mm	
Coaxial Focusing Device	Coarse focusing moving range:25mm Fine focusing moving range: 2μm/division	
Abbe Condenser	N.A. 1.25 Abbe condenser	
Diaphragm	Iris diaphragm	
Illumination	3W LED Aspheric illumination	
Input Voltage	Wide 100V-240V voltage	
Plug	European standard, American standard, Germany Standard ,Australian standard, British standard etc.	

Product Application



About Us



Chongqing Scope Instrument Co., Ltd. is a professional supplier who focus on laboratory instrument, and our team has more than 10 years rich experience in the field of laboratory instrument. We have perfect relationship with most of leading manufacturers in the field of laboratory instrument in China, and supply perfect service and different kinds of quality laboratory instrument from different leading manufacturers to our customers. We now serving for more than 80 customers from Europe, Asia, Africa, South America and Australia. Also, we have established good relationship with more than 200 customers in worldwide market. Our laboratory instrument are widely used in the field of labs, education, medicine, petroleum, geological mineral and chemical industry, and electronic industry.

Product Application



Payment

Western Union

T/T



Shipment



Are you a factory or trading company??

We are trading company. But we have established long-term relationships with factories. And we provide high quality products for you.

Can we get support if we have our own market position?

Please inform us your detailed mind on your market demand. We will discuss and propose helpful suggestions for you, to find the best solution for you.

If I place an order, how long it will take before I got delivery??

It would largely depend on the size of order and the items you want.

Can I use our own logo??

Yes, customized your logo is available.

Do you provide warranty for the goods?

Yes, we provide 1 year warranty for all the goods from us.

